



*A Trusted Supplier of
Ophthalmic OPD & O.T. EYE Equipments*

I Care Industry

i Care Industry is a leading trading firm of Diagnostic Optical Instruments, medical furniture and hospital equipments. Our major customers consist of Government Hospitals, Private Hospitals, domestic and overseas dealers.

Our wide range of products consist of Auto Referecatometer, Lensometer, Vision Drums, Acuity Charts, Trail Lens Sets, Operation Tables, Refraction Chair Unit, Zoom Microscopes, Different Eye Testing Equipments and other miscellaneous products.

We strictly follow Fair Business Practices, believes in transparency and ethics. We are committed to providing our customers Economical, Cost Effective, High Quality Products and Services in the shortest possible time schedules.

We are committed to preserving our natural resources, we strictly monitor our processes to minimize the impact on our environment.

Quality Assurance

We maintain strict quality parameters at all levels of our operations and manufacturing process such as recruitment and retaining technically competent persons, procurement of quality raw-materials and components, Technically Advanced Machineries, R & D, prompt after-sales services, timely delivery and execution of orders. We are ISO 9001-2015 Certified / CE Compliant company.


PRODUCTS GROUP	PAGE NO.
SPECULUMS	1 - 3
RETRACTORS * KNIVES	4
KERATOMES * KNIVES * SPUDS	5 - 6
CURETTES * PROBES * SPOONS * DIALATORS	7 - 8
DISSECTORS * CLAMPS	9
MALLETS * FORCEPS	10
NIBBLER * PUNCHES * TREPHINES	11
HOOKS	12 - 13
CAUTERYS * LOOPS * SPATULAS	14 - 15
SPATULAS * CALIPERS * MARKERS	16
MARKERS * FIXATION RING	17
FORCEPS	18 - 33
HOLDER & FOLDER	34
I.O.L. MANIPULATORS	35 - 36
HOOKS * POLISHERS * NUCLEUS ITEMS	37
CHOPPERS	38
VITREO-RETINAL INSTRUMENT	39
SCISSORS	40 - 45
BLADE BREAKERS * NEEDLE HOLDERS	46 - 48
CANNULAE * HANDLE	49 - 56
TREPHINE PUNCH	57
STERILIZATION CASE	57
BIMANUAL SETS & HANDLE	58

CLEANING

- Distilled or mineral water (Lab tested soft water) should be used to Decontaminate & rinse the instruments.
- For removal of contaminated stuck on the surface, use soft bristled brush or ultrasonic cleaner.
- Cannulas and hand pieces should be cleaned with distilled water.
- Then use syringe or compressed dry air to remove wet.

STERILIZATION

- Use distilled water in autoclave.
- Autoclave should be rustproof from inside otherwise rust may transfer by steam.
- During dry heat sterilization the temperature should be fully controlled otherwise instruments hardness may get affected.



BARRAQUER

Wire speculum

ME-0006 (A) adult

ME-0006 (B) child


ME-0006 (C) large



BARRAQUER

Wire eye speculum


ME-0012 medium



GRANDON-BARRAQUER

Wire speculum


ME-0018



KRATZ-BARRAQUER

Wire speculum

ME-0024



LIEBERMAN

Eye speculum

ME-0030 (A)

Phaco Speculum

ME-0030 (A)



CLARK

Eye speculum

ME-0036 (A) adult

ME-0036 (B) child



LANCASTER

Eye speculum

ME-0042



MURDOCK

Eye speculum

ME-0048 (A) Small

ME-0048 (B) Medium


ME-0048 (C) Large



SAUER

Eye speculum

ME-0054



INFANT

Wire speculum

ME-0060




WEISS

Eye speculum

ME-0066 (A) Adult

ME-0066 (B) Child



WILLIAMS

Eye speculum

ME-0072 (A) adult

ME-0072 (B) child



SCLERAL DEPRESSOR

ME-0078 (A) single

ME-0078 (B) double



SCHEPENS

Scleral depressor

ME-0084 (A) 8mm


ME-0084 (B) 10mm



MUELLER'S

Lacrimal Sac Retractor

ME-0090



STEVENSON

Lacrimal sac retractor

ME-0096



JAFFE

Lid retractor

ME-0102 (A) small

ME-0102 (B) large



DESMARRES
Lid retractor
ME-0108



PAUL
Retractor, blunt
ME-0114



CATPAW
Retractor 4
prongs for sac
ME-0120



SCHEPENS
Orbital retractor
ME-0126



DASTOOR
Iris retractor
ME-0132



DASTOOR
Iris retractor
ME-0138 right & left



DASTOOR
Iris retractor
ME-0144



BARRAQUER

Cataract knife in
sliding case

ME-0156



SCHNAUDI GL

Hockey knife

ME-0162



GRAEFE

Cataract knife

ME-0150



BARRAQUER

Keratoplasty Knife

ME-0198



ZIEGLER

Iris knife

ME-0174 (A) 3mm

ME-0174 (B) 5mm

ME-0174 (C) 6mm



TOOK KNIFE

Corneal knife

ME-0180



GILLS

Corneal knife

ME-0186



PAUFIQUE

Keratoplasty Knife

ME-0192

**BARRAQUER**

Keratoplasty Knife
ME-0198

**BOWMAN**

Discission needle
ME-0204 (A) straight
ME-0204 (B) curved

**MEYER**

Schwickerath localizer
ME-0210

**GRAEFE**

Cystotme
ME-0216

**GRAEFE**

Cystotme & spoon
ME-0222

**DIX**

Foreign body spud
ME-0228

**ROLLET**

Rugine for lacrimal sac
ME-0234

**GRAEFE**

Scarifier for lacrimal sac
ME-0240



MEYHOEFER

Chalazion curette sharp

ME-0246 (A) 1mm

ME-0246 (B) 2mm

ME-0246 (C) 3mm



HARMS

Trabeculotomy probe


ME-0252 Right & Left



WILDER

Lacrimal dilator

ME-0258



BOWMAN

Lacrimal Probe set

ME-0264 (4 pes set)



MULE

Evisceration spoon

ME-0270



WELLS

Enucleation spoon

ME-0276



BUNGE
Evisceration
spoon
ME-0282



GRAEFE
Lens spoon
ME-0288



WEST
Bone Chisel
ME-0294



WEST
Bone gouge
ME-0300



FREER
Twin elevator
sharp & blunt
ME-0306



TRAQUIAR'S

Periostoral elevator
double-end

ME-0312



DASTOOR

Lacrimal sac dissector

ME-0318



LANG

Lacrimal sac dissector

ME-0324



PIGTAIL PROBE

With suture hole

ME-0330



SERREFINE

ME-0336 (A) straight

ME-0336 (B) curved



BACKHAUS

Towel clamp


ME-0342



TOWEL CLAMP

Cross action

ME-0348



MALLET
For DCR
ME-0354



HARTMANN
Mosquito forcep 3.5"
delicate
ME-0360 (A) Straight
ME-0360 (B) Curved



HALSTED
Mosquito hemostatic
forceps 5.5"
ME-0366 (A) Straight
ME-0366 (B) Curved



MALLET

For DCR
ME-0354



HARTMANN

Mosquito forcep 3.5"
delicate

ME-0360 (A) Straight
ME-0360 (B) Curved



HALSTED

Mosquito hemostatic
forceps 5.5"

ME-0366 (A) Straight
ME-0366 (B) Curved



KERRISON

D C R Punch

ME-0372 (A) 1.0mm

ME-0372 (B) 2.0mm

ME-0372 (C) 3.0mm

ME-0372 (D) 4.0mm



KELLY'S

Detachment Punch

ME-0378



GASS

Trabeculotomy punch

ME-0384



CASTROVIEJO

Disposable trephine

ME-0390 5mm to 9mm



GRAEFE

Iris hook

ME-0396 (A) sharp

ME-0396 (B) blunt



ARRUGA

Hook

ME-0402



CULLER

Muscle hook

ME-0408



GREEN

Strabismus hook

ME-0414




JAMESON

Strabismus hook

ME-0420 (A) small

ME-0420 (B) large



GRAEFE

Strabismus hook

ME-0426



O'CONNOR
Tenotomy hook
ME-0432



GASS
Detachment hook
With oval hole
ME-0438



FLAT
Hook
ME-0444



CHAVASSE'S
Muscle hook
ME-0450



SMITH
Lens expressor
ME-0456



KIRBY
Lens expressor
Hook double
ME-0462



LENS EXPRESSOR

Hook & Loop

ME-0468



KIRBY

Lens expressor
hook & spatula

ME-0474



WILLS

Hospital cautery
cooper ball

ME-0480



WILLS

Cautery ball point ss

ME-0486



JOSHI

Cautery copper Tip


ME-0492



BLUNT

Cautery

ME-0498



DR. PARUL SHAH'S

Type bimanual set

ME-2202

Simcoe Cortex 22 Gauge
Teplon Handle Knob
with Tubing



DR. J.K. SHAH'S

Type bimanual set

ME-2208

Teplon Handle
Simcoe Cortex
1/A (23g) wash
Cannula (19g)



I/A SYSTEM

Irrigating aspirating
ss with titanium tip

ME-2214



BIMANUAL

Handpiece for bimanual
technique (T)

ME-2220 (A) Irrigating
ME-2220 (B) Aspirating



B.P. HANDLE

ME-2226



TREPHINE PUNCH
ME-2172




S S BOX

Metal tray with one silicone mat
ME-2178

Metal tray with two silicone mat
ME-2184

Metal tray with four silicone mat
ME-2190



P C BOX

ME-2196
Small 6x2.5x0.75
Medium 7.5x4x0.75
Large 10x6x0.75



KELMAN

Cystotome with sharp point and sharp cutting edge 22gauge
ME-2136



KNOLLE-KELMAN

Cystotome small blade with sharp point and sharp cutting edge 22gauge
ME-2142



Mc INTYRE

Cystotome angled 10mm from tip with sharp point, 23gauge
ME-2148



SILICONE BULB

With adaptor
ME-2154



TEPLON HANDLE

ME-2160



INFUSION

Handle

ME-2166 (A) 9.0 Long
ME-2166 (B) 12.5 Long



SIMCOE

Cortex extractor cannula
smooth closed tip, with side
opening facing into curved
23 gauge sandblasted

ME-2100 (A) 13mm

ME-2100 (B) 21mm



SIMCOE

Cortex extractor and capsule
polisher front opening facing
into curved 23gauge sandblasted

ME-2106 (A) 13mm

ME-2106 (B) 21mm



ANIS

Cortex aspirating cannula with
0.4mm side port finely
sandblasted tip 23gauge

ME-2112 (A) Straight

ME-2112 (B) Angled



GRAETHER

Collar button Iris retracting &
capsular polisher cannula

ME-2118 (A) Straight

ME-2118 (B) Angled



DALJIT

Blast cannula

ME-2124



HESLIN

Infusion cannula with silicone
tubing and adaptor

ME-2130 (A) 23 Gauge

ME-2130 (B) 26 Gauge



DREWS

Posterior capsule polisher slightly angled, 2mm i.d. loop with sandblasted base 0.2 irrigating port facing into circle 25gauge
[ME-2064](#)



JENSEN

Posterior capsule polisher carbide impregnated on posterior surface of olive tip front opening gently curved
[ME-2070 \(A\) 25 Gauge](#)
[ME-2070 \(B\) 23 Gauge](#)



JENSEN

Capsule polisher carbide impregnated all around olive tips front opening gently curved
[ME-2076 \(A\) 25 Gauge](#)
[ME-2076 \(B\) 23 Gauge](#)



KNOLLE

Posterior capsule polisher 2mm angled tip, front opening finely sandblasted
[ME-2082 \(A\) 23 Gauge](#)
[ME-2082 \(B\) 30 Gauge](#)



KRATZ

Scratcher front opening with carbide impregnated tips slightly angled 22gauge
[ME-2088](#)



SHAPIRO

Posterior capsule polisher carbide impregnated tip with dual 0.3mm side ports curved 23gauge
[ME-2094](#)



KNOLLE

Anterior chamber irrigating cannula 2mm angled tip front opening

ME-2022 (A) 30 Gauge
ME-2022 (B) 23 Gauge



GIMBEL

Hydro-dissection Cannula smooth blunt tip front opening 30gauge

ME-2028



AIR INJECTION

Cannula angled 5mm from tip

ME-2034 (A) 26 Gauge
ME-2034 (B) 27 Gauge
ME-2034 (C) 30 Gauge



TROUTMAN

Alpha-chymotrypsin cannula delicate olive tip 25gauge

ME-2040



KUGLEN

Irrigating lens manipulator disc shaped tip with 0.25mm diameter peg protruding from its center, 25gauge

ME-2046 (A) Straight
ME-2046 (B) Angled



YOUENS

Lens manipulating cannula front opening angled 30gauge

ME-2052



WELSH

Flat olive tip cannula 25gauge

ME-2058 (A) Straight
ME-2058 (B) Right
ME-2058 (C) Left



GILLS

Irrigating-aspirating
cannula side by side
front opening with
silicone tubing 23gauge
ME-1980



'J' SHAPED

Irrigating aspirating Cannula
ME-1986 (A) 12 'O' Clock
ME-1986 (B) Right
ME-1986 (C) Left



McINTYRE

Inner cannula 23gauge
ME-1992



McINTYRE

Outer cannula straight
blunt tip front
opening 19gauge
ME-1998



McINTYRE

Anterior chamber cannula
smooth blunt tip
ME-2004 23 Gauge Angled
ME-2010 26 Gauge Angled
ME-2016 30 Gauge Angled



KNOLLE-PEARCE

Irrigating vectis three front opening ports
23 gauge loop


ME-1944



DREWS

Irrigating vectis
one front and two backward opening parts, 23 gauge loop

ME-1950



SIMCOE

Irrigating lens
nucleus loop two front opening ports finely, serrated, strong curved 25 gauge

ME-1956 (A) Right

ME-1956 (B) Left



SHEETS

Irrigating vectis three large front opening ports
serrated 20 gauge loop

ME-1962



SIMCOE

Irrigating aspirating
cannula 23 gauge
with silicone tubing

ME-1968 Regular



SIMCOE

Irrigating-aspirating
cannula with silicone
tubing 23 gauge

ME-1974 Reverse



MORRISON

Aspirating cannula
smooth blunt tip
angled at 7mm 23gauge
ME-1902



VISCOELASTIC

Aspirating cannula
angled at 10mm from tip
ME-1908 (A) 19 Gauge
ME-1908 (B) 22 Gauge



BRACKEN

Irrigating cannula
smooth beveled tip
front opening 19 gauge
ME-1914



BINKHORST

Irrigating cannula
smooth blunt tip
front opening 22 gauge
MS-1920 (A) straight
MS-1920 (B) right
MS-1920 (C) left



WILBRANDT-BINKHORST

Irrigating cannula smooth
blunt tip front
opening 23 gauge
ME-1926 (A) Straight
ME-1926 (B) Right
ME-1926 (C) Left



MCINTRYE- BINKHORST

Irrigating cannula smooth
blunt tip front
opening 26 gauge
ME-1932 (A) Straight
ME-1932 (B) Right
ME-1932 (C) Left



WEST

Irrigating cannula
blunt tip, 23 gauge
side opening
ME-1938



LACRIMAL

Cannula, 23gauge
ME-1836 (A) Straight
ME-1836 (B) Curved



CANNULA

A.C. Wash curved 20gauge
ME-1842



BISHOP-HARMON

A.C. Wash cannula 10mm
19 gauge
ME-1848



RANDOLPH

Cyclodialysis cannula
curved spatulated tip
angled at 13mm,
ME-1854 (A) 19 Gauge
ME-1854 (B) 21 Gauge



CASTROVIEJO

Cyclodialysis cannula
top opening
ME-1860 (A) 10mm
ME-1860 (B) 15mm




KELMAN

Cyclodialysis cannula
spatulated tip side
opening 19gauge
ME-1866 (A) Right
ME-1866 (B) Left



SALZ

Aspirating cannula
15mm long 23gauge
straight thin wall
ME-1872



GILLS

Aspirating cannula
smooth blunt tip
angled 25 gauge
ME-1878 5mm Tip
ME-1884 7mm Tip
ME-1890 8mm Tip
ME-1896 10mm Tip



BARRAQUER

Micro needle holder medium

ME-1794 (A) without Lock

ME-1794 (B) with Lock



BARRAQUER

Micro needle holder short model

ME-1800 (A) without Lock

ME-1800 (B) with Lock



BARRAQUER

Micro needle holder extra delicate

ME-1806 (A) Straight without Lock

ME-1806 (B) Straight with Lock

ME-1806 (C) Curved without Lock

ME-1806 (D) Curved with Lock



NEEDLE HOLDER

Micro extra delicate

ME-1812 (A) Straight without Lock

ME-1812 (B) Straight with Lock

ME-1812 (C) Curved without Lock

ME-1812 (D) Curved with Lock



TROUTMAN

Micro needle holder

ME-1818 (A) without Lock

ME-1818 (B) with Lock



MCINTYRE

Needle holder

ME-1824 (A) without Lock

ME-1824 (B) with Lock



LIMS

Needle holder

ME-1830 (A) without Lock

ME-1830 (B) with Lock



CASTROVIEJO

Needle holder standard jaws

ME-1752 (A) Straight without Lock

ME-1752 (B) Straight with Lock

ME-1752 (C) Curved without Lock

ME-1752 (D) Curved with Lock



CASTROVIEJO

Micro needle holder delicate jaws

ME-1758 (A) Straight without Lock

ME-1758 (B) Straight with Lock

ME-1758 (C) Curved without Lock

ME-1758 (D) Curved with Lock



JAFFE

Needle holder delicate jaws

ME-1764 (A) Straight without Lock

ME-1764 (B) Straight with Lock

ME-1764 (C) Curved without Lock

ME-1764 (D) Curved with Lock



Mc PHERSON

Needle holder delicate jaws

ME-1770 (A) Straight without Lock

ME-1770 (B) Straight with Lock

ME-1770 (C) Curved without Lock

ME-1770 (D) Curved with Lock



BARRAQUER

Needle holder standard model

ME-1776 (A) without Lock

ME-1776 (B) with Lock



BARRAQUER

Needle holder

ME-1782 (A) without Lock

ME-1782 (B) with Lock




BARRAQUER

Needle holder

ME-1788 (A) without Lock

ME-1788 (B) with Lock



SWISS

Model blade holder
pencil 12 cm.

ME-1722




SWISS

Model blade holder
pencil type

ME-1728 (A) Straight 9 cm

ME-1728 (B) Angled 9 cm



BARRAQUER

Blade breaker & holder

ME-1734




KALT

Needle holder

ME-1740 (A) Straight

ME-1740 (B) Curved



SILCOCK

Needle holder

MS-1746




CASTROVIEJO

Corneal scissor
blunt tips large blade
ME-1686 Slight Cur ved



CASTROVIEJO

Corneal scissor
blunt tips large blade
ME-1692 Half Cur ved



CASTROVIEJO

Corneal scissor
blunt tips large blade
ME-1698 Full Curved



CASTROVIEJO

Corneal scissor
ME-1704 (A) Right Small
ME-1704 (B) Left Small



ME-1710 (A) Right Medium
ME-1710 (B) Left Medium

ME-1716 (A) Right Large
ME-1716 (B) Left Large



MICRO

Corneal scissor

blunt tips small blade

ME-1656 Full Curved



CASTROVIEJO

Corneal scissor

blunt tips large blade

ME-1662 (A) Slightly Curved

ME-1662 (B) Half Curved

ME-1662 (C) Full Curved



MICRO

Corneal scissor

blunt tips small blade (hole)

ME-1668 Slightly Curved



MICRO

Corneal scissor

blunt tips small blade (hole)

ME-1674 Half Curved



MICRO

Corneal scissor

blunt tips small blade (hole)

ME-1680 Full Curved



WESTCOTT
Stitch scissor
standard blade
ME-1620



WESTCOTT
Stitch scissor
sharp pointed tips
medium blade
ME-1626



CASTROVIEJO
Corneal scissor
angled to side blunt tips
ME-1632 (A) Small Blade 11mm
ME-1632 (B) Large Blade 15mm



IRIS
Scissor sharp tips
small blade
ME-1638 (A) Straight
ME-1638 (B) Curved



MICRO
Corneal scissor
blunt tips small blade
ME-1644 Slight Cur ved



MICRO
Corneal Scissor
blunt tips small blade
ME-1650 Half Cur ved



GILLS - VANNAS

Capsulotomy scissor sharp tips

ME-1584

Angled Forward 6mm & 11mm



STERN-GILLS

Capsulotomy scissor sharp tips


ME-1590 Curved



STERN-GILLS

Capsulotomy scissor
pointed blade

ME-1596 Angled 45° 10mm Tips



WESTCOTT

Tenotomy scissor curved
blunt tips standard blade

ME-1602 (A) Right

ME-1602 (B) Left



WESTCOTT

Tenotomy scissor curved
blunt tips medium blade


ME-1608



WESTCOTT

Conjunctival scissor curved
blunt tips small blade

ME-1614



STEVENS

Tenotomy scissor

ME-1548 (A) Straight

ME-1548 (B) Curved




NOYES

Iris scissor

ME-1554 (A) Blunt

ME-1554 (B) Sharp



VANNAS

Scissor sharp tips

ME-1560 (A) Straight Blade 7mm

ME-1560 (B) Straight Blade 10mm



VANNAS

Scissor sharp tips

ME-1566 (A) Curved Blade 7mm

ME-1566 (B) Curved Blade 10mm



VANNAS

Scissor sharp tips

ME-1572 Angled to Side



CAPSULOTOMY

Scissor curved blunt tips

ME-1578 14mm Blade



ENUCLEATION

Scissor

ME-1524 (A) Half Cur ved

ME-1524 (B) Full Cur ved

ME-1524 (C) ST Blunt Tip



IRIS

Scissor 9 cm

ME-1530 (A) Straight

ME-1530 (B) Cur ved



IRIS

Scissor 10.5 cm

ME-1536 (A) Straight

ME-1536 (B) Cur ved




STRABISMUS

ME-1542 (A) Straight

ME-1542 (B) Cur ved


DORC TYPE FORCEPS & SCISSOR



ME-1410
Straight Forceps



ME-1416
Serrated Forceps



ME-1422
End Gripping



ME-1428
45° Forceps



ME-1434
Forceps



ME-1440
Straight Scissor
Sharp Tip



ME-1446
Curved Scissor
Sharp Tip



ME-1452
Vertical Scissor
45°



ME-1458
Vertical Scissor
90°



ME-1464
Capsularhexis
45°




ME-1470
Capsularhexis
90°



ME-1476
Foreign Body

MEMBRANE PEELAR




ME-1482
Peeler



ME-1488
Spatula Curved
Pointed Tip



ME-1494
Scratcher Claw
Type



ME-1500
Spatula Hooked



ME-1506
Spatula Hooked
Short Length



ME-1512
Scratcher



ME-1518
Spatula Knob



KANSAS

Nucleus vectis, solid
used as a "cutting block"

ME-1368 (A) Right
ME-1368 (B) Left



GILLS-WELSH

Lens loop with olive tip

ME-1374



ATRAUMATIC

Ring forcep with
sinskey hook

ME-1380



PHACO CHOPPER

Universal

ME-1386



PHACO CHOPPER

Ball tip

ME-1392



PHACO CHOPPER

Nagahara

ME-1398



PRECHOPPER WITH

Substainer

ME-1404 (A) Straight

ME-1404 (B) Angled

ME-1404 (C) Straight (Hole)

ME-1404 (D) Angled (Hole)



HIRSCHMAN

Iris push & pull hook
ME-1320 (A) Straight
ME-1320 (B) Angled



RENTSCH "BOAT HOOK"

Lens manipulator
ME-1326



TWIST HOOK

For scleral fixation
ME-1332 (A) Right
ME-1332 (B) Left



NIGHTINGALE

Capsule polishing
curette 1.7mm
ME-1338



NIGHTINGALE

Capsule polishing Curette
ME-1344 (A) Right
ME-1344 (B) Left



NORDAN

Capsule elevator
for in-th-bag
IOL placement
ME-1350



KANSAS

Nucleus trisector
ME-1356 (A) Right
ME-1356 (B) Left



KANSAS

Nucleus bisector
ME-1362



KUGLEN

Iris hook and lens manipulator
push-pull model

ME-1266 (A) Straight

ME-1266 (B) Angled



FENZL

Lens manipulating hook

T shaped tip

ME-1272 (A) Straight

ME-1272 (B) Angled



LESTER

Lens manipulator

ME-1278 (A) Straight

ME-1278 (B) Angled



KNOLLE

Iris hook very blunt

ME-1284 Straight



BONN

Micro iris hook


ME-1290



SHEPARD

Iris hook U-shaped

ME-1296



SIEPSEr

Endocapsular controller

trumpet shaped tip

ME-1302



PISACANO

Nucleus rotator

flat blunt tip

ME-1308



MANN-HUNKELER

Capsule retractor

with 0.3mm diameter
semi-ball

ME-1314



BECHERT

Lens nucleus rotator
blunt tip angled
ME-1206 (A) 7mm from Tip
ME-1206 (B) 9mm from Tip



CLAYMAN

Nucleus rotator
v-notched tip
ME-1212



HIRSCHMAN

Lens spatula
with small notches
ME-1218



SHEETS

Lens spatula
flat smooth
ME-1224



JAFFE

Lens spatula
flat smooth
ME-1230



GRAETHER

Collar button
micro iris retractor
ME-1236 (A) Straight
MS-1236 (B) Angled



JAFFE-KNOLLE

Iris hook very blunt
small angled
ME-1242



SINSKEY

Lens manipulating hook
ME-1248 (A) Straight
ME-1248 (B) Angled



LEWICKY

Lens manipulating hook
angled
ME-1254



DRYSDALE

Nucleus manipulator
paddle shaped tip
ME-1260



LENS
Injector
ME-1158



LENS
Injector
ME-1164



SILICON
Forcep
ME-1170 (A) Folder
ME-1170 (B) Inserter




LIVERNOIS
Silicone lens forcep
ME-1176 Folder



MCDONALD
I.O.L Inserting forcep
ME-1182




FOLDABLE
Forcep
ME-1188 (A) Folder
ME-1188 (B) Inserter



FAULKNER
Forcep
ME-1194 (A) Folder
ME-1194 (B) Inserter



BURRATTO
Soft I.O.L Implantation
ME-1200



KNOLLE-SHEPARD

Lens holding forcep
with peg serrated jaws
ME-1128



SHEPARD

Lens holding forcep
with serrated lower jaws
ME-1134 (A) with Lock
ME-1134 (B) without Lock



SHEPARD II

Lens holding forcep
with nudge notch
ME-1140 (A) with Lock
ME-1140 (B) without Lock



CLAYMAN

Lens holding forcep
delicate angled jaws
ME-1146 (A) with Lock
ME-1146 (B) without Lock




MOORE

Lens introducing forcep
ME-1152 (A) with Lock
ME-1152 (B) without Lock



COLIBRI
Ring tip forcep
ME-1092



RING TIP FORCEP
Long model angled
ME-1098



ARRUGA
Capsular forcep
ME-1104 Standard



MICRO
Ring forcep
ME-1110



BLAYDES
Lens holding forcep
ME-1116



ANIS
Lens holding forcep
extra delicate
ME-1122



BORES

Incision spreading forcep
ME-1062



WATZKE

Sleeve spreading forcep
ME-1068



DODICK

Nucleus cracker
cross-action
ME-1074



KANSAS

Nucleus removal forcep
ME-1080



DALJIT

Lens holding forcep
ME-1086



KELMAN-Mc PHERSON

Forcep

ME-1032 (A) 10mm Straight Tying

ME-1032 (B) 10mm Angled Tying




KELMAN-Mc PHERSON

Forceps 11mm plain

Angled delicate

ME-1038



KRAFF

Capsule tag forcep

11mm long jaws

ME-1044



UTRATA

Capsularhexis forcep

delicate grasping tips

ME-1050



MASKET

Capsularhexis forcep

ME-1056



Mc PHERSON

Forcep

ME-1002 (A) 8mm Straight Tying

ME-1002 (B) 8mm Angled Tying



Mc PHERSON

Forcep

ME-1008 (A) 10mm Straight Tying

ME-1008 (B) 10mm Angled Tying



Mc PHERSON

Forcep delicate 11mm angled

ME-1014



KELMAN-Mc PHERSON

Forcep

ME-1020 (A) 5mm Straight Tying

ME-1020 (B) 5mm Angled Tying



KELMAN-Mc PHERSON

Forcep

ME-1026 (A) 8mm Straight Tying

ME-1026 (B) 8mm Angled Tying



HEATH

Chalazion forcep
ME-0972 12x14mm



FRANCIS

Chalazion forcep
ME-0978 12x14mm Oval



DESMARRES

Chalazion forcep
ME-0984 (A) 20mm
ME-0984 (B) 27mm
ME-0984 (C) 31mm



SWISS

Model muscle forcep
ME-0990 (A) Right
ME-0990 (B) Left



Mc PHERSON

Forcep
ME-0996 (A) 5mm Straight Tying
ME-0996 (B) 5mm Angled Tying



SNELLEN

Entropion forcep

ME-0942 (A) Right Small
ME-0942 (B) Left Small
ME-0942 (C) Right Large
ME-0942 (D) Left Large



BERKE

Ptosis clamp

ME-0948 (A) 20mm
ME-0948 (B) 27mm



AYER

Chalazion forcep

ME-0954



LAMBERT

Chalazion forcep


ME-0960



BAIRD

Chalazion forcep

ME-0966 8x12mm




BONN

Micro tying forcep

ME-0906 (A) Short

ME-0906 (B) Long



JAYLE'S

Tying forcep

ME-0912 Straight



JEWELERS

Forcep

ME-0918



JEWELERS

Forcep angled

ME-0924 (A) Straight

ME-0924 (B) Curved



SATTLER

Advancement forcep

ME-0930 (A) Right

ME-0930 (B) Left




WORTH

Advancement forcep

ME-0936 (A) Right

ME-0936 (B) Left



HARMS

Tying forcep smooth jaws

ME-0870 (A) Straight

ME-0870 (B) Curved



SINSKEY

Tying forcep curved

ME-0876



TENNANT

Tying forcep extra delicate jaws

ME-0882 (A) Straight

ME-0882 (B) Curved




MICRO

Micro tying forcep

ME-0888 (A) Straight Small

ME-0888 (B) Curved Small



MICRO

Micro tying forcep

ME-0894 (A) Straight Large

ME-0894 (B) Straight Large


ME-0894 (C) Angled



CASTROVIEJO

Wide handle tying forcep

ME-0900



PIERSE

Corneal forcep angled

ME-0834



ELSCHMIG

Capsular forcep

ME-0840



SWISS MODEL

Extra capsular forcep

ME-0846



BARRAQUER

Cilia forcep slightly angled


ME-0852



ZIEGLER

Cilia forcep

ME-0858



BARRAQUER

Tying forcep


ME-0864 (A) Straight
ME-0864 (B) Curved



TROUTMAN

Superior-rectus forcep

ME-0798



PIERSE TYPE

Micro forcep

ME-0804 (A) 1.5mm Tip

ME-0804 (B) 0.2mm Tip




PIERSE TYPE

Micro forcep curved

ME-0810(A) 0.2mm Tip

ME-0810 (B) 0.5mm Tip



PIERSE TYPE

Fine forcep

ME-0816 (A) 0.2mm Tip

ME-0816 (B) 0.3mm Tip

ME-0816 (C) 0.8mm Tip



PIERSE

Colibri Micro forcep

ME-0822 0.2mm Tip



PIERSE

Corneal forcep delicate

ME-0828



GRAEFE

Fixation forcep

ME-0762



BARRAQUER

Fixation forcep

ME-0768



BORES

Corneal fixation
forcep 3mm 1x2 teeth

ME-0774 (A) Straight

ME-0774 (B) Angled



ARROWSMITH

Corneal fixation
forcep 8 mm 1x2 teeth

ME-0780 (A) Straight Lock

ME-0780 (B) Straight without Lock

ME-0780 (C) Angled with Lock

ME-0780 (D) Angled without Lock



KRAMER

Fixation forcep 13mm1x2 teeth

ME-0786 (A) Straight Lock

ME-0786 (B) Straight without Lock

ME-0786 (C) Angled with Lock


ME-0786 (D) Angled without Lock



DASTOOR

Superior-rectus forcep


ME-0792



CASTROVIEJO

Wide handle 1x2
teeth forcep

ME-0726



PAUFIQUE

Suture forcep 1x2 teeth

ME-0732



ALABAMA

Suture forcep

ME-0738



JAYLE'S

Suture forcep 1x2 teeth


ME-0744



MOORFIELD

Suture and conjunctival
forcep

ME-0750



FIXATION

Forcep

ME-0756 (A) 1x2 Teeth

ME-0756 (B) 2x3 Teeth

ME-0756 (C) 3x4 Teeth



CASTROVIEJO

Colibri forcep
very delicate 1x2 teeth

ME-0690



TENNANT

Colibri corneal forcep

ME-0696



GIRARD

Colibri forcep 1x2 teeth
with tying platform

ME-0702



LIM'S

Corneal forcep 1x2 teeth

ME-0708 (A) Small

ME-0708 (B) Large



HOFMANN-POLACK

Corneal suturing forcep
with 1x2 pierce teeth

ME-0714



CASTROVIEJO


Suture forcep teeth

ME-0720 (A) 0.12 mm Straight

ME-0720 (B) 0.5 mm Straight

ME-0720 (C) Angled

ME-0720 (D) 0.5mm ST with Lock



BISHOP-HARMON

Iris forcep

ME-0654 (A) Extra Delicate

ME-0654 (B) 1x2 Teeth

ME-0654 (C) Serrated



KIRBY

Iris forcep 1x2 teeth


ME-0660



UTILITY

Forcep


ME-0666 1x2 Straight



BONACCOL-TO

Utility forcep


ME-0672



WILLS HOSPITAL

Utility forcep

ME-0678



BARRAQUER

Colibri forcep

ME-0684 (A) 0.12mm Teeth


ME-0684 (B) 0.2mm Teeth



EYE DRESSING

Forcep serrated


ME-0618 (A) Straight
ME-0618 (B) Curved



IRIS FORCEP

1x2 teeth

ME-0624 (A) Straight
ME-0624 (B) Curved



SWISS

Model iris forcep
curved 1x2 teeth


ME-0630



HESS

Iris forcep


ME-0636 (A) Swiss Model
ME-0636 (B) Standard Model



BONN

Iris forcep

ME-0642 (A) with Platform
ME-0642 (B) without Platform



BONN

Iris forcep


ME-0648 (A) with Platform
ME-0648 (B) without Platform



HOFFER

Optic zone marker
With cross hairs

ME-0582 (A) 3.0mm
ME-0582 (B) 3.5mm
ME-0582 (C) 4.0mm
ME-0582 (D) 4.5mm
ME-0582 (E) 5.0mm



THORNTON

Optic zone marker
With cross pin

ME-0588 (A) 3.0mm
ME-0588 (B) 3.5mm
ME-0588 (C) 4.0mm
ME-0588 (D) 4.5mm
ME-0588 (E) 5.0mm



GRENA

Visual axis marker

ME-0594



MARKER

ME-0600 (A) Double
ME-0600 (B) Triple




CASTROVIEJO

Caliper 20 mm

ME-0564 (A) Straight


ME-0564 (B) Curved



STAHL

Caliper

ME-0570



OSHER-NEUMANN

Radial marker

ME-0576 (A) 4 Blades

ME-0576 (B) 6 Blades

ME-0576 (C) 8 Blades

ME-0576 (D) 12 Blades

ME-0576 (E) 16 Blades



THORNTON

Fixation ring with

Atraumatic teeth 16 mm

ME-0606



HOFMAN-THORNTON

Global fixation ring with

ME-0612 (A) 16mm

ME-0612 (B) 13mm




SPATULA
With needle
ME-0546



PATON
Keratoplasty spatula
ME-0552 sharp & blunt



CASTROVIEJO
Cyclodialysis spatula
ME-0558 (A) Double
ME-0558 (B) Elsching single
ME-0558 (C) Spatula with
retractor



SMITH-FISHER
Iris spatula double
ME-0528



WECKER
Iris spatula
ME-0534



GREEN
Spatula double end
ME-0540

CATARACT SET

1.	0006	WIRE SPECULUM	01
2.	0258	LACRIMAL DIALATOR	07
3.	0264	PROBE SET	07
4.	0336	SERREFINE (2)	09
5.	0342	BACKHAUS TOWEL CLAMP (2)	09
6.	0360	HARTMANN MOSQUITO FORCEP (2)	10
7.	0444	FLAT HOOK	13
8.	0456	SMITH LENS EXPRESSOR	13
9.	0480	WILLS HOSPITAL CAUTERY COPPER BALL	14
10.	0504	SNELLEN LENS LOOPS	15
11.	0516 B	IRIS REPOSITOR DOUBLE	15
12.	0564 A	CASTROVIEJO CALIPER	16
13.	0648 A	BONN IRIS FORCEP	18
14.	0684	BARRAQUER COLIBRI FORCEP	19
15.	0792	DASTOOR SUPERIOR FORCEP	22
16.	0870 B	HARMS TYING FORCEP	25
17.	1002 B	Mc PHERSON FORCEP 8mm. ANG	29
18.	1026 B	KELMAN MC PHERSON FORCEP 8mm. ANG	29
19.	1050	UTRATA CAPSULAR HEXIS FORCEP	30
20.	1248 B	SINSKEY LENS MANIPULATING HOOK	35
21.	1566 B	VANNAS SCISSOR	41
22.	1608	WESTCOTT TENOTOMY SCISSOR	42
23.	1650	MICRO CORNEAL SCISSOR HALF CUV	43
24.	1740 A	KALT NEEDLE HOLDER	46
25.	1794 A	BARRAQUER NEEDLE HOLDER W/O LOCK	48
26.	1848	ANTERIOR CHAMBER WASH CANNULA 27g	49
27.	1968	SIMCOE IRRIGATING-ASPIRATING CANNULA	51
28.	2034 B	AIR INJECTION CANNULA	53
29.	2154	SILICONE BULB WITH ADAPTOR	56

IOL SET

1.	0006	WIRE SPECULUM	01
2.	0336	SERREFINE (2)	09
3.	0342	BACKHAUS TOWEL CLAMP (2)	09
4.	1008 B	MC PHERSON FORCEP 10mm. ANG	29
5.	1032 B	KELMAN MC PHERSON FORCEP 10mm. ANG	30
6.	1044	RSFF CAPSULE TAG FORCEP	30
7.	1050	UTRATA CAPSULAR HEXSIS FORCEP	30
8.	1116	BLADES LENS HOLDING FORCEP	32
9.	1236 B	GRAETHER COLLAR BUTTON	35
10.	1248 B	SINSKEY LENS MANIPULATOR HOOK	35
11.	1266 B	KUGLEN IRIS HOOK	36
12.	1338	NIGHTINGALE CAPSULE POLISHER	37
13.	1968	SIMCOE IRRIGATING-ASPIRATING CANNULA	51
14.	2034 B	AIR INJECTION CANNULA	53
15.	2154	SILICONE BULB WITH ADAPTOR	56
16.	2166 A	IRRIGATING-ASPIRATING HANDLE	56

CHALAZION SET

1.	0006	WIRE SPECULUM	01
2.	0246 B	MEYHOEFER CHALAZION CURETTE	07
3.	0342	BACKHAUS TOWEL CLAMP (2)	09
4.	0360 A	HARTMANN MOSQUITO FORCEP (2)	10
5.	0954	AYER CHALAZION FORCEP	27
6.	0960	LAMBERT CHALAZION FORCEP	27
7.	0984 B	DESMARRES CHALAZION FORCEP	28
8.	1722	BLADE HOLDER PENCIL MODEL	46
9.	2154	SILICONE BULB WITH ADAPTOR	56
10.	2226	BARD PARKER HANDLE	58

ENUCLEATION SET

1.	0042	LANCASTER EYE SPECULUM	02
2.	0276	WELLS ENUCLEATION SPOON	07
3.	0360 B	HARTMANN MOSQUITO FORCEP	10
4.	0396 A	GRAEFE MUSCLE HOOK	12
5.	0756 A	FIXATION FORCEP 1x2	21
6.	1524 A	ENUCLEATION SCISSOR	40
7.	1740 A	KALT NEEDLE HOLDER	46

SQUINT SET

1.	0036 A	CLARK SPECULUM	02
2.	0336 A	SERREFINE (2)	09
3.	0444	FLAT HOOK	13
4.	0564 A	CASTROVIEJO CALIPER	16
5.	0936	ADVANCEMENT FORCEP RIGHT & LEFT	26
6.	1542 A	STRABISMUS SCISSOR	40
7.	1548 B	TENOTOMY SCISSOR	41
8.	1602 A	WESTCOTT SCISSOR	42

GLAUCOMA SURGERY

1.	0006	WIRE SPECULUM	01
2.	0252	HARMS TRABECULOTOME PROBE RIGHT & LEFT	07
3.	0342	BACKHAUS TOWEL CLAMP (2)	09
4.	0360	HARTMANN MOSQUITO FORCEP (2)	10
5.	0384	GASS CORNEOSCLERAL PUNCH	11
6.	0558 A	CASTROVIEJO CYCLODIALYSIS SPATULA D/E	16
7.	0564 A	CASTROVIEJO CALIPER	16
8.	0618 A	EYE DRESSING FORCEP SERRATED	18
9.	0654 A	BISHOP-HARMON IRIS FORCEP	19
10.	0684	BARRAQUER COLIBRI FORCEP	19
11.	0870 A	HARMS TYING FORCEP	25
12.	1032 B	KELMAN MC PHERSON FORCEP 10mm. ANG	30
13.	1530 A	IRIS SCISSOR	40
14.	1566	VANNAS SCISSOR	41
15.	1602 A	WESTCOTT TENOTOMY SCISSOR	42

16.	1620	WESTCOTT STITCH SCISSOR	43
17.	1752 C	CASTROVIEJO NEEDLE HOLDER W/O LOCK	47
18.	1794 A	BARRAQUER NEEDLE HOLDER W/O LOCK	48
19.	2034 B	AIR INJECTION CANNULA 27g	53
20.	2154	SILICONE BULB WITH ADAPTOR	56

LID SET

1.	0108	LID RETRACTOR	04
2.	0246 B	CHALAZION CURETTE	07
3.	0336 B	SERREFINE (2)	09
4.	0342	BACKHAUS TOWEL CLAMP (2)	09
5.	0360 A	HARTMANN MOSQUITO FORCEP (4)	10
6.	0426	GRAEFE MUSCLE HOOK	12
7.	0654 B	BISHOP-HARMON FORCEP	18
8.	0762	GRAEFE FIXATION FORCEP	22
9.	0870 A	HARMS TYING FORCEP	25
10.	0948 A	PTOSIS CLAMP	27
11.	1026 B	KELMAN Mc PHERSON FORCEP 8mm. ANG	29
12.	1530 A	IRIS SCISSOR	40
13.	1548 A	TENOTOMY SCISSOR	41
14.	1620	WESTCOTT TENOTOMY SCISSOR	43
15.	1626	WESTCOTT STITCH SCISSOR	43
16.	1794 A	BARRAQUER NEEDLE HOLDER W/O LOCK	48
17.	2154	SILICONE BULB WITH ADAPTOR	56

SAC & DCR SET

1.	0096	STEVENSON LACRIMAL SAC RETRACTOR	03
2.	0120	CATPAW RETRACTOR (2)	04
3.	0258	LACRIMAL DIALATOR	07
4.	0264	LACRIMAL PROBE SET	07
5.	0294	BONE CHISEL	08
6.	0300	BONE GOUGE	08
7.	0312	TRAQUAIR'S PERIOSTERAL ELEVATOR	09
8.	0318	DASTOOR LACRIMAL SAC DETECTOR	09
9.	0360 A	HARTMANN MOSQUITO FORCEP (4)	10
10.	0372	BONE PUNCH	11
11.	1740 A	KALT NEEDLE HOLDER	46
12.	2154	SILICONE BULB WITH ADAPTOR	56

VITREO RETINAL SET

1.	0036 A	CLARK SPECULUM	02
2.	0084 A	SCHEPENS SCLERAL DEPRESSOR	03
3.	0126	SCHEPENS ORBITAL RETRACTOR	04
4.	0336 A	SERREFINE	09
5.	0342	BACKHAUS TOWEL CLAMP (2)	09
6.	0360 A	HARTMANN MOSQUITO FORCEP (4)	10
7.	0420 A	JAMESON STRABISMUS HOOK	12
8.	0438	GASS DETACHMENT HOOK	13
9.	0564 A	CASTROVIEJO CALIPER	16
10.	0666 B	UTILITY FORCEP SERRATED	19

11.	0684	BARRAQUER COLIBRI FORCEP	19
12.	0732	PAUFIQUE SUTURE FORCEP	21
13.	0870	HARMS TYING FORCEP	25
14.	1068	WATZKE SPREADING FORCEP	31
15.	1416	V R SERRATED FORCEP	39
16.	1422	V R END GRIPPING	39
17.	1446	V R CURVED SCISSOR	39
18.	1452	VERTICAL 45O SCISSOR	39
19.	1482	MEMBRANE PEELAR	39
20.	1566 B	VANNAS IRIS SCISSOR	41
21.	1608	WESTCOTT TENOTOMY SCISSOR	42
22.	1758 C	CASTROVIEJO NEEDLE HOLDER CUV W/O LOCK	47

KERATOPLASTY SET

1.	0006	WIRE SPECULUM	01
2.	0366	ARTRY FORCEP (ST/CR) (6)	10
3.	0390	TREPHINE	11
4.	0516 A	IRIS REPOSITOR	15
5.	0684	COLIBRI FORCEP	19
6.	0708	LIMS FORCEP	20
7.	0816	HOSKIN FORCEP	23
8.	0870 A	SUTURE TYING FORCEP S/T	25
9.	0870 B	SUTURE TYING FORCEP C/R	25
10.	1530	RING SCISSOR (ST/CR)	40
11.	1584	GILLS VANNAS SCISSOR	42
12.	1608	WESTCOTT SCISSOR	42
13.	1674	CORNEAL SCISSOR	44
14.	1704	LEFT & RIGHT SCISSOR	45
15.	1830	NEEDLE HOLDER	48
16.	1896	A.C.WASH CANNULA	49
17.	1968	SIMCOE CANNULA	51
18.	2172	TREPHINE PUNCH	57